

For Research Use Only

REST Recombinant antibody, PBS Only (Capture)

Catalog Number: 84846-6-PBS



Basic Information

Catalog Number:

84846-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17707

GenBank Accession Number:

BC132859

GeneID (NCBI):

5978

UNIPROT ID:

Q13127

Full Name:

RE1-silencing transcription factor

Calculated MW:

1097 aa, 122 kDa

Purification Method:

Protein A purification

CloneNo.:

241789G1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01629-3, Human REST Recombinant Matched Antibody Pair - PBS only. 84846-6-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag17707. 84846-2-PBS was HRP conjugated as the detection antibody. Range: 19.5-1250 pg/mL



The mean REST concentration was determined to be 801.98 ng/mL in Raji cell extract based on a 1.0 mg/mL extract load.