

For Research Use Only

FOXP4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84855-3-PBS



Basic Information

Catalog Number:

84855-3-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10047

GenBank Accession Number:

BC052803

GeneID (NCBI):

116113

UNIPROT ID:

Q8IVH2

Full Name:

forkhead box P4

Calculated MW:

678 aa, 73 kDa

Purification Method:

Protein A purification

CloneNo.:

242356E3

Applications

Tested Applications:

FC (Intra), Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01621-2, Human FOXP4 Recombinant Matched Antibody Pair - PBS only. 84855-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag10047. 84855-3-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean FOXP4 concentration was determined to be 5.4 ng/mL in HepG2 cell extract based on a 1.0 mg/mL extract load.