

For Research Use Only

CHRNA3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84860-2-PBS



Basic Information

Catalog Number:

84860-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0272

GenBank Accession Number:

BC000513

GeneID (NCBI):

1136

UNIPROT ID:

P32297

Full Name:

cholinergic receptor, nicotinic, alpha 3

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

242410A8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01633-2, CHRNA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84860-2-PBS. Detection antibody: 84860-3-PBS. Standard: Ag0272. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01633-1, CHRNA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84860-2-PBS. Detection antibody: 84860-1-PBS. Standard: Ag0272. Range: 0.781-100 ng/mL.