

For Research Use Only

# RIG-1/DDX58 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84861-1-PBS



## Basic Information

Catalog Number:

84861-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18585

GenBank Accession Number:

BC132786

GeneID (NCBI):

23586

UNIPROT ID:

O95786

Full Name:

DEAD (Asp-Glu-Ala-Asp) box  
polypeptide 58

Calculated MW:

925 aa, 106 kDa

Purification Method:

Protein A purification

CloneNo.:

242424A10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

# Selected Validation Data



Cytometric bead array standard curve of MP01632-1, RIG-1/DDX58 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84861-1-PBS. Detection antibody: 84861-3-PBS. Standard: Ag18585. Range: 0.625-80 ng/mL.