

For Research Use Only

RIG-1/DDX58 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84861-4-PBS



Basic Information

Catalog Number:

84861-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18585

GenBank Accession Number:

BC132786

GeneID (NCBI):

23586

UNIPROT ID:

O95786

Full Name:

DEAD (Asp-Glu-Ala-Asp) box
polypeptide 58

Calculated MW:

925 aa, 106 kDa

Purification Method:

Protein A purification

CloneNo.:

242424G10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01632-2, RIG-1/DDX58 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84861-4-PBS. Detection antibody: 84861-2-PBS. Standard: Ag18585. Range: 0.625-80 ng/mL.