

For Research Use Only

# CENPF Recombinant antibody, PBS Only (Capture)

Catalog Number: 84865-3-PBS



## Basic Information

Catalog Number:

84865-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG29665

GenBank Accession Number:

NM\_016343

GeneID (NCBI):

1063

UNIPROT ID:

P49454

Full Name:

centromere protein F, 350/400ka  
(mitosin)

Calculated MW:

368 kDa

Purification Method:

Protein A purification

CloneNo.:

242277H7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01637-2, CENPF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84865-3-PBS. Detection antibody: 84865-1-PBS. Standard: Ag29665. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01637-1, CENPF Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84865-3-PBS. Detection antibody: 84865-2-PBS. Standard: Ag29665. Range: 0.625-80 ng/mL