

For Research Use Only

VWA5A Recombinant antibody, PBS Only (Capture)

Catalog Number: 84868-3-PBS



Basic Information

Catalog Number:

84868-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22198

GenBank Accession Number:

BC001234

GeneID (NCBI):

4013

UNIPROT ID:

O00534

Full Name:

von Willebrand factor A domain containing 5A

Calculated MW:

786 aa, 86 kDa

Purification Method:

Protein A purification

CloneNo.:

242320G8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01639-1, VWA5A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84868-3-PBS. Detection antibody: 84868-1-PBS. Standard: Ag22198. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01639-2, VWA5A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84868-3-PBS. Detection antibody: 84868-2-PBS. Standard: Ag22198. Range: 0.781-100 ng/mL.