

For Research Use Only

# ARRDC1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84874-4-PBS



## Basic Information

Catalog Number:

84874-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15341

GenBank Accession Number:

BC032346

GeneID (NCBI):

92714

UNIPROT ID:

Q8N5I2

Full Name:

arrestin domain containing 1

Calculated MW:

433 aa, 46 kDa

Purification Method:

Protein A purification

CloneNo.:

242396E7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01649-1, ARRDC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84874-4-PBS. Detection antibody: 84874-3-PBS. Standard: Ag15341. Range: 0.625-80 ng/mL.