

For Research Use Only

NNMT Recombinant antibody, PBS Only (Capture)

Catalog Number: 84876-4-PBS



Basic Information

Catalog Number:

84876-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7197

GenBank Accession Number:

BC000234

GeneID (NCBI):

4837

UNIPROT ID:

P40261

Full Name:

nicotinamide N-methyltransferase

Calculated MW:

30 kDa

Purification Method:

Protein A purification

CloneNo.:

242394E4

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01628-3, Human NNMT Recombinant Matched Antibody Pair - PBS only. 84876-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag7197. 84876-3-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL



The mean NNMT concentration was determined to be 94.15 ng/mL in L02 cell extract based on a 1.30 mg/mL extract load and 139.70 ng/mL in HeLa cell extract based on a 1.00 mg/mL extract load.