

For Research Use Only

NDUFS2 Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84880-3-PBS



Basic Information

Catalog Number:

84880-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28008

GenBank Accession Number:

BC000170

GeneID (NCBI):

4720

UNIPROT ID:

O75306

Full Name:

NADH dehydrogenase (ubiquinone)

Fe-S protein 2, 49kDa (NADH-coenzyme Q reductase)

Calculated MW:

53 kDa

Purification Method:

Protein A purification

CloneNo.:

242238H7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01653-1, NDUF52 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84880-2-PBS. Detection antibody: 84880-3-PBS. Standard: Ag28008. Range: 0.938-60 ng/mL.



Cytometric bead array standard curve of MP01653-2, NDUF52 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84880-3-PBS. Detection antibody: 84880-1-PBS. Standard: Ag28008. Range: 0.938-60 ng/mL.