

For Research Use Only

CHKA Recombinant antibody, PBS Only (Capture)

Catalog Number: 84881-1-PBS



Basic Information

Catalog Number:

84881-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4761

GenBank Accession Number:

BC036471

GeneID (NCBI):

1119

UNIPROT ID:

P35790

Full Name:

choline kinase alpha

Calculated MW:

439 aa, 52 kDa

Purification Method:

Protein A purification

CloneNo.:

242095B11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01644-2, CHKA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84881-1-PBS. Detection antibody: 84881-2-PBS. Standard: Ag4761. Range: 1.563-200 ng/mL.



Cytometric bead array standard curve of MP01644-1, CHKA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84881-1-PBS. Detection antibody: 84881-3-PBS. Standard: Ag4761. Range: 1.563-200 ng/mL.