

For Research Use Only

GRIA3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84885-1-PBS



Basic Information

Catalog Number:

84885-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30898

GenBank Accession Number:

NM_007325

GeneID (NCBI):

2892

UNIPROT ID:

P42263

Full Name:

glutamate receptor, ionotropic,
AMPA 3

Calculated MW:

101 kDa

Purification Method:

Protein A purification

CloneNo.:

242392B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01656-1, Glutamate receptor 3/GluA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84885-1-PBS. Detection antibody: 84885-2-PBS. Standard: Ag30898. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01656-2, Glutamate receptor 3/GluA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84885-1-PBS. Detection antibody: 84885-3-PBS. Standard: Ag30898. Range: 0.781-100 ng/mL