

For Research Use Only

Mouse Galectin-3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84889-1-PBS



Basic Information

Catalog Number: 84889-1-PBS	GenBank Accession Number: GeneID (NCBI): 16854	Purification Method: Protein A purification
Size: 1 mg/ml	UNIPROT ID: P16110	CloneNo.: 242453C5
Source: Rabbit	Full Name: lectin, galactose binding, soluble 3	
Isotype: IgG	Calculated MW: 28kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01641-2, MOUSE Galectin-3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84889-2-PBS. Detection antibody: 84889-1-PBS. Standard: Eg1232. Range: 0.469-60 ng/mL.