

For Research Use Only

AMACR Recombinant antibody, PBS Only (Detector)

Catalog Number: 84891-6-PBS



Basic Information

Catalog Number:

84891-6-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4005

GenBank Accession Number:

NM_014324.6

GeneID (NCBI):

23600

UNIPROT ID:

Q9UHK6-1

Full Name:

alpha-methylacyl-CoA racemase

Calculated MW:

42 kDa

Purification Method:

Protein A purification

CloneNo.:

250503E1

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01623-4, Human AMACR/p504S Recombinant Matched Antibody Pair - PBS only. 84891-7-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg4005. 84891-6-PBS was HRP conjugated as the detection antibody. Range: 62.5-4000 pg/mL



The mean AMACR/p504S concentration was determined to be 2,611.4 pg/mL in LNCaP cell extract based on a 2.6 mg/mL extract load and 3,344.8 pg/mL in MCF-7 cell extract based on a 1.3 mg/mL extract load.