

For Research Use Only

# PMS2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84894-2-PBS



## Basic Information

Catalog Number:

84894-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33508

GenBank Accession Number:

BC093921

GeneID (NCBI):

5395

UNIPROT ID:

P54278

Full Name:

PMS2 postmeiotic segregation increased 2 (S. cerevisiae)

Calculated MW:

862 aa, 96 kDa

Purification Method:

Protein A purification

CloneNo.:

242455B12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01663-1, PMS2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84894-2-PBS. Detection antibody: 84894-1-PBS. Standard: Ag33508. Range: 0.781-100 ng/mL.