

For Research Use Only

# TRIM56 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84896-1-PBS



## Basic Information

Catalog Number:

84896-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG21992

GenBank Accession Number:

BC048194

GeneID (NCBI):

81844

UNIPROT ID:

Q9BRZ2

Full Name:

tripartite motif-containing 56

Calculated MW:

755 aa, 81 kDa

Purification Method:

Protein A purification

CloneNo.:

242497A7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01657-1, TRIM56 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84896-1-PBS. Detection antibody: 84896-2-PBS. Standard: Ag21992. Range: 0.781-100 ng/mL.