## For Research Use Only

## TRIM56 Recombinant antibody

Catalog Number:84896-7-RR



**Basic Information** 

Catalog Number: 84896-7-RR Concentration:

1000 µ g/ml
Source:
Rabbit
Isotype:

IgG tripartite motif-containing 56

GenBank Accession Number:

BC048194

81844

Q9BRZ2 Full Name:

81 kDa

GeneID (NCBI):

UNIPROT ID:

Immunogen Catalog Number:Calculated MW:AG21992755 aa, 81 kDaObserved MW:

Tested Applications:

Species Specificity:

human

**ELISA** 

Applications

## **Background Information**

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer


PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

Purification Method: Protein A purfication

CloneNo.: 242497D3

## Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed MCF-7 cells using TRIM56 antibody (84896-7-RR, Clone: 242497D3) at dilution of 1:200 and Coralite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).