

For Research Use Only

TCOF1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84914-1-PBS



Basic Information

Catalog Number:

84914-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1447

GenBank Accession Number:

BC014559

GeneID (NCBI):

6949

UNIPROT ID:

Q13428

Full Name:

Treacher Collins-Franceschetti
syndrome 1

Calculated MW:

97 kDa, 144 kDa, 152 kDa

Purification Method:

Protein A purification

CloneNo.:

242201D8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01672-1, TCOF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84914-3-PBS. Detection antibody: 84914-1-PBS. Standard: Ag1447. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01672-2, TCOF1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84914-2-PBS. Detection antibody: 84914-1-PBS. Standard: Ag1447. Range: 0.625-80 ng/mL