

For Research Use Only

FEN1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84916-4-PBS



Basic Information

Catalog Number:

84916-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6552

GenBank Accession Number:

BC000323

GeneID (NCBI):

2237

UNIPROT ID:

P39748

Full Name:

flap structure-specific endonuclease 1

Calculated MW:

43 kDa

Purification Method:

Protein A purification

CloneNo.:

242330B3

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01694-3, Human FEN1 Recombinant Matched Antibody Pair - PBS only. 84916-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag6552. 84916-2-PBS was HRP conjugated as the detection antibody. Range: 0.781-50 ng/mL



The mean FEN1 concentration was determined to be 112.14 ng/mL in Raji cell extract based on a 1.00 mg/mL extract load and 70.29 ng/mL in HepG2 cell extract based on a 1.00 mg/mL extract load.