

For Research Use Only

SHOC2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84917-4-PBS



Basic Information

Catalog Number:

84917-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG11754

GenBank Accession Number:

BC050445

GeneID (NCBI):

8036

UNIPROT ID:

Q9UQ13

Full Name:

soc-2 suppressor of clear homolog (C. elegans)

Calculated MW:

582 aa, 65 kDa

Purification Method:

Protein A purification

CloneNo.:

242343F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01674-2, SHOC2/Sur-8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84917-2-PBS. Detection antibody: 84917-4-PBS. Standard: Ag11754. Range: 0.781-100 ng/mL.