

For Research Use Only

HECTD1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84918-4-PBS



Basic Information

Catalog Number:

84918-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14581

GenBank Accession Number:

BC006237

GeneID (NCBI):

25831

UNIPROT ID:

Q9ULT8

Full Name:

HECT domain containing 1

Calculated MW:

2612 aa, 290 kDa

Purification Method:

Protein A purification

CloneNo.:

242344H7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01675-2, HECTD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84918-2-PBS. Detection antibody: 84918-4-PBS. Standard: Ag14581. Range: 0.781-100 ng/mL.