

For Research Use Only

Mouse CD25/IL-2RA Recombinant antibody, PBS Only (Detector)

Catalog Number: 84925-2-PBS



Basic Information

Catalog Number:

84925-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_008367.3

GeneID (NCBI):

16184

UNIPROT ID:

P01590

Full Name:

interleukin 2 receptor, alpha chain

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

242494F8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01669-2, MOUSE CD25/IL2RA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84925-1-PBS. Detection antibody: 84925-2-PBS. Standard: Eg3087. Range: 0.625-80 ng/mL