

For Research Use Only

# PRICKLE1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84926-3-PBS



## Basic Information

Catalog Number:

84926-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18407

GenBank Accession Number:

BC114939

GeneID (NCBI):

144165

UNIPROT ID:

Q96MT3

Full Name:

prickle homolog 1 (Drosophila)

Calculated MW:

831 aa, 94 kDa

Purification Method:

Protein A purification

CloneNo.:

242506H4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01670-2, PRICKLE1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84926-3-PBS. Detection antibody: 84926-1-PBS. Standard: Ag18407. Range: 0.156-20 ng/mL



Cytometric bead array standard curve of MP01670-1, PRICKLE1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84926-3-PBS. Detection antibody: 84926-2-PBS. Standard: Ag18407. Range: 0.156-20 ng/mL