

For Research Use Only

# TRIM27 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84927-1-PBS



## Basic Information

Catalog Number: 84927-1-PBS	GenBank Accession Number: BC013580	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 5987	CloneNo.: 242538A11
Source: Rabbit	UNIPROT ID: P14373	
Isotype: IgG	Full Name: tripartite motif-containing 27	
Immunogen Catalog Number: AG2849	Calculated MW: 513 aa, 58 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01671-2, TRIM27 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84927-1-PBS. Detection antibody: 84927-2-PBS. Standard: Ag2849. Range: 0.313-40 ng/mL