

For Research Use Only

TRIM27 Recombinant antibody, PBS Only

Catalog Number: 84927-5-PBS



Basic Information

Catalog Number:

84927-5-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2849

GenBank Accession Number:

BC013580

GeneID (NCBI):

5987

UNIPROT ID:

P14373

Full Name:

tripartite motif-containing 27

Calculated MW:

513 aa, 58 kDa

Observed MW:

58 kDa

Purification Method:

Protein A purification

CloneNo.:

242538B5

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84927-5-RR (TRIM27 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84927-5-PBS in a different storage buffer formulation.