

For Research Use Only

# TEAD4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84932-5-PBS



## Basic Information

Catalog Number: 84932-5-PBS	GenBank Accession Number: BC015497	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 7004	CloneNo.: 242496F4
Source: Rabbit	UNIPROT ID: Q15561	
Isotype: IgG	Full Name: TEA domain family member 4	
Immunogen Catalog Number: AG3059	Calculated MW: 434 aa, 48 kDa	

## Applications

Tested Applications:  
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP01682-3, Human TEAD4 Recombinant Matched Antibody Pair - PBS only. 84932-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag3059. 84932-4-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL



The mean TEAD4 concentration was determined to be 8.65 ng/mL in HepG2 cell extract based on a 1.00 mg/mL extract load and 6.31 ng/mL in MCF-7 cell extract based on a 1.20 mg/mL extract load.