

For Research Use Only

# RPS3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84947-1-PBS



## Basic Information

Catalog Number:

84947-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7346

GenBank Accession Number:

BC003137

GeneID (NCBI):

6188

UNIPROT ID:

P23396

Full Name:

ribosomal protein S3

Calculated MW:

27 kDa

Purification Method:

Protein A purification

CloneNo.:

242504A9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01686-1, RPS3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84947-1-PBS. Detection antibody: 84947-2-PBS. Standard: Ag7346. Range: 0.313-40 ng/mL.