

For Research Use Only

LARP7 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84950-1-PBS



Basic Information

Catalog Number:

84950-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10850

GenBank Accession Number:

BC006981

GeneID (NCBI):

51574

UNIPROT ID:

Q4G0J3

Full Name:

La ribonucleoprotein domain family, member 7

Calculated MW:

582 aa, 67 kDa

Purification Method:

Protein A purification

CloneNo.:

242561C5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01697-2, LARP7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84950-3-PBS. Detection antibody: 84950-1-PBS. Standard: Ag10850. Range: 1.563-100 ng/mL.