

For Research Use Only

# GOLGA3 Recombinant antibody

Catalog Number: 84955-4-RR



## Basic Information

Catalog Number:

84955-4-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15531

GenBank Accession Number:

BC142658

GeneID (NCBI):

2802

UNIPROT ID:

Q08378

Full Name:

golgi autoantigen, golgin subfamily a, 3

Calculated MW:

1498 aa, 167 kDa

Observed MW:

167 kDa

Purification Method:

Protein A purification

CloneNo.:

242603C7

Recommended Dilutions:

WB 1:1000-1:8000

## Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : HEK-293 cells, HeLa cells, L02 cells, PC-3 cells

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84955-4-RR (GOLGA3 antibody) at dilution of 1:4000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 84955-4-RR against Human GOLGA3 were performed. The affinity constant is 0.475 nM.