

For Research Use Only

# RCN2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 84956-1-PBS



## Basic Information

Catalog Number:

84956-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0247

GenBank Accession Number:

BC004892

GeneID (NCBI):

5955

UNIPROT ID:

Q14257

Full Name:

reticulocalbin 2, EF-hand calcium binding domain

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

242612B3

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01713-2, RCN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84956-1-PBS. Detection antibody: 84956-2-PBS. Standard: Ag0247. Range: 0.469-60 ng/mL