

For Research Use Only

# TFB1M Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 84957-1-PBS



## Basic Information

Catalog Number:

84957-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9785

GenBank Accession Number:

BC017788

GeneID (NCBI):

51106

UNIPROT ID:

Q8WVM0

Full Name:

transcription factor B1, mitochondrial

Calculated MW:

346 aa, 40 kDa

Purification Method:

Protein A purification

CloneNo.:

242524A5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01712-2, TFB1M Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84957-2-PBS. Detection antibody: 84957-1-PBS. Standard: Ag9785. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01712-1, TFB1M Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84957-1-PBS. Detection antibody: 84957-3-PBS. Standard: Ag9785. Range: 0.781-100 ng/mL.