

For Research Use Only

# SUGT1 Recombinant antibody, PBS Only

Catalog Number: 84963-5-PBS



## Basic Information

Catalog Number:

84963-5-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2275

GenBank Accession Number:

BC000911

GeneID (NCBI):

10910

UNIPROT ID:

Q9Y2Z0

Full Name:

SGT1, suppressor of G2 allele of SKP1  
(S. cerevisiae)

Calculated MW:

41 kDa

Observed MW:

38-41 kDa

Purification Method:

Protein A purification

CloneNo.:

242590G3

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84963-5-RR (SUGT1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84963-5-PBS in a different storage buffer formulation.