

For Research Use Only

HNRNPU Recombinant antibody, PBS Only (Capture)

Catalog Number: 84964-1-PBS



Basic Information

Catalog Number:

84964-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6141

GenBank Accession Number:

BC003621

GeneID (NCBI):

3192

UNIPROT ID:

Q00839

Full Name:

heterogeneous nuclear
ribonucleoprotein U (scaffold
attachment factor A)

Calculated MW:

825 aa, 91 kDa

Purification Method:

Protein A purification

CloneNo.:

242609B4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01690-2, HNRNPU Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84964-1-PBS. Detection antibody: 84964-3-PBS. Standard: Ag6141. Range: 0.469-60 ng/mL.