

For Research Use Only

PBK/SPK Recombinant antibody, PBS Only (Detector)

Catalog Number: 84970-2-PBS



Basic Information

Catalog Number:

84970-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9033

GenBank Accession Number:

BC015191

GeneID (NCBI):

55872

UNIPROT ID:

Q96KB5

Full Name:

PDZ binding kinase

Calculated MW:

322 aa, 36 kDa

Purification Method:

Protein A purification

CloneNo.:

242545C2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01693-2, PBK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84970-3-PBS. Detection antibody: 84970-2-PBS. Standard: Ag9033. Range: 0.781-100 ng/mL.