

For Research Use Only

C3orf20 Recombinant antibody, PBS Only

Catalog Number: 84975-1-PBS



Basic Information

Catalog Number:

84975-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG15031

GenBank Accession Number:

BC038406

GeneID (NCBI):

84077

UNIPROT ID:

Q8ND61

Full Name:

chromosome 3 open reading frame 20

Calculated MW:

904 aa, 101 kDa

Observed MW:

88 kDa

Purification Method:

Protein A purification

CloneNo.:

242404E2

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

The function of C3orf20 has not been widely studied, and is yet to be fully elucidated. Two isoforms were produced by alternative splicing with predicted MW of 88 and 101 kDa according to UniProt. Catalog#20906-1-AP can recognize both 88 and 101 kDa isoforms.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 84975-1-RR (C3orf20 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 84975-1-PBS in a different storage buffer formulation.