

For Research Use Only

# YBX1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 84993-1-PBS



## Basic Information

Catalog Number:

84993-1-PBS

Concentration:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14129

GenBank Accession Number:

BC010430

GeneID (NCBI):

4904

UNIPROT ID:

P67809

Full Name:

Y box binding protein 1

Calculated MW:

324 aa, 36 kDa

Purification Method:

Protein A purification

CloneNo.:

242462D10

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01766-2, YBX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 84993-2-PBS. Detection antibody: 84993-1-PBS. Standard: Ag14129. Range: 1.563-200 ng/mL.