

For Research Use Only

GRK5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85015-1-PBS



Basic Information

Catalog Number:

85015-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10511

GenBank Accession Number:

BC064506

GeneID (NCBI):

2869

UNIPROT ID:

P34947

Full Name:

G protein-coupled receptor kinase 5

Calculated MW:

590 aa, 68 kDa

Purification Method:

Protein A purification

CloneNo.:

242490A7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01747-2, GRK5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85015-2-PBS. Detection antibody: 85015-1-PBS. Standard: Ag10511. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01747-1, GRK5 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85015-3-PBS. Detection antibody: 85015-1-PBS. Standard: Ag10511. Range: 0.781-100 ng/mL