

For Research Use Only

NDST1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85016-2-PBS



Basic Information

Catalog Number:

85016-2-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG24422

GenBank Accession Number:

BC012888

GeneID (NCBI):

3340

UNIPROT ID:

P52848

Full Name:

N-deacetylase/N-sulfotransferase
(heparan glucosaminyl) 1

Calculated MW:

882 aa, 101 kDa

Purification Method:

Protein A purification

CloneNo.:

242509D10

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01768-2, NDST1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85016-2-PBS. Detection antibody: 85016-3-PBS. Standard: Ag24422. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01768-1, NDST1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85016-2-PBS. Detection antibody: 85016-1-PBS. Standard: Ag24422. Range: 0.313-40 ng/mL