

For Research Use Only

# DOCK11 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85022-1-PBS



## Basic Information

Catalog Number:

85022-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG35831

GenBank Accession Number:

NM\_144658

GeneID (NCBI):

139818

UNIPROT ID:

Q5JSL3

Full Name:

dedicator of cytokinesis 11

Calculated MW:

238 kDa

Purification Method:

Protein A purification

CloneNo.:

242575C9

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01759-1, DOCK11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85022-2-PBS. Detection antibody: 85022-1-PBS. Standard: Ag35831. Range: 0.781-100 ng/mL



Cytometric bead array standard curve of MP01759-2, DOCK11 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85022-3-PBS. Detection antibody: 85022-1-PBS. Standard: Ag35831. Range: 0.781-100 ng/mL