

For Research Use Only

FXYD5 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85025-1-PBS



Basic Information

Catalog Number:

85025-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2808

GenBank Accession Number:

BC009642

GeneID (NCBI):

53827

UNIPROT ID:

Q96DB9

Full Name:

FXYD domain containing ion transport
regulator 5

Calculated MW:

178 aa, 19 kDa

Purification Method:

Protein A purification

CloneNo.:

242581A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01770-2, Dysadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85025-2-PBS. Detection antibody: 85025-1-PBS. Standard: Ag2808. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01770-1, Dysadherin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85025-3-PBS. Detection antibody: 85025-1-PBS. Standard: Ag2808. Range: 0.313-40 ng/mL