

For Research Use Only

PDLIM1/CLP36 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85027-1-PBS



Basic Information

Catalog Number:

85027-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2274

GenBank Accession Number:

BC000915

GeneID (NCBI):

9124

UNIPROT ID:

O00151

Full Name:

PDZ and LIM domain 1

Calculated MW:

329 aa, 36 kDa

Purification Method:

Protein A purification

CloneNo.:

242585C8

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01749-1, PDLIM1/CLP36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85027-3-PBS. Detection antibody: 85027-1-PBS. Standard: Ag2274. Range: 0.313-40 ng/mL



Cytometric bead array standard curve of MP01749-2, PDLIM1/CLP36 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85027-2-PBS. Detection antibody: 85027-1-PBS. Standard: Ag2274. Range: 0.313-40 ng/mL