

For Research Use Only

# HEXIM1 Recombinant antibody, PBS Only

Catalog Number: 85028-3-PBS



## Basic Information

Catalog Number:

85028-3-PBS

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8144

GenBank Accession Number:

BC006460

GeneID (NCBI):

10614

UNIPROT ID:

O94992

Full Name:

hexamethylene bis-acetamide  
inducible 1

Calculated MW:

41 kDa

Purification Method:

Protein A purification

CloneNo.:

242587F1

## Applications

Tested Applications:

Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 85028-3-PBS against Human HEXIM1 were performed. The affinity constant is 0.181 nM.