

For Research Use Only

HEXIM1 Recombinant antibody, PBS Only (Detector)

Catalog Number:85028-4-PBS



Basic Information

Catalog Number: 85028-4-PBS	GenBank Accession Number: BC006460	Purification Method: Protein A purification
Concentration: 1000 ug/ml	GeneID (NCBI): 10614	CloneNo.: 242587D5
Source: Rabbit	UNIPROT ID: O94992	
Isotype: IgG	Full Name: hexamethylene bis-acetamide inducible 1	
Immunogen Catalog Number: AG8144	Calculated MW: 41 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP01740-3, Human HEXIM1 Recombinant Matched Antibody Pair - PBS only. 85028-5-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag8144. 85028-4-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL



The mean HEXIM1 concentration was determined to be 122.36 ng/mL in HeLa cell extract based on a 1.50 mg/mL extract load, 285.55 ng/mL in MCF-7 cell extract based on a 3.50 mg/mL extract load and 310.04 ng/mL in Jurkat cell extract based on a 3.10 mg/mL extract load.