

For Research Use Only

C1QC Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85031-1-PBS



Basic Information

Catalog Number:

85031-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG10400

GenBank Accession Number:

BC009016

GeneID (NCBI):

714

UNIPROT ID:

P02747

Full Name:

complement component 1, q
subcomponent, C chain

Calculated MW:

245 aa, 26 kDa

Purification Method:

Protein A purification

CloneNo.:

242624A2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01746-1, C1qC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85031-1-PBS. Detection antibody: 85031-2-PBS. Standard: Ag10400. Range: 0.781-100 ng/mL.



Cytometric bead array standard curve of MP01746-2, C1qC Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85031-3-PBS. Detection antibody: 85031-1-PBS. Standard: Ag10400. Range: 0.781-100 ng/mL.