

For Research Use Only

# RFK Recombinant antibody, PBS Only (Capture/Detector)

Catalog Number: 85032-3-PBS



## Basic Information

Catalog Number:

85032-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG8591

GenBank Accession Number:

BC007069

GeneID (NCBI):

55312

UNIPROT ID:

Q969G6

Full Name:

riboflavin kinase

Calculated MW:

162 aa, 18 kDa

Purification Method:

Protein A purification

CloneNo.:

242628G5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01741-1, RFK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85032-2-PBS. Detection antibody: 85032-3-PBS. Standard: Ag8591. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01741-2, RFK Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85032-3-PBS. Detection antibody: 85032-1-PBS. Standard: Ag8591. Range: 0.625-80 ng/mL