

For Research Use Only

Mouse Slamf7 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85033-1-PBS



Basic Information

Catalog Number:

85033-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2672

GenBank Accession Number:

NM_144539.5

GeneID (NCBI):

75345

UNIPROT ID:

Q8BHK6-1

Full Name:

SLAM family member 7

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

242629B5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01769-2, MOUSE Slanf7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85033-1-PBS. Detection antibody: 85033-3-PBS. Standard: Eg2672. Range: 0.625-80 ng/mL



Cytometric bead array standard curve of MP01769-1, MOUSE Slanf7 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85033-1-PBS. Detection antibody: 85033-2-PBS. Standard: Eg2672. Range: 0.625-80 ng/mL