

For Research Use Only

# SEC22B Recombinant antibody, PBS Only (Capture)

Catalog Number: 85037-3-PBS



## Basic Information

Catalog Number: 85037-3-PBS	GenBank Accession Number: BC001364	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 9554	CloneNo.: 242659E4
Source: Rabbit	UNIPROT ID: O75396	
Isotype: IgG	Full Name: SEC22 vesicle trafficking protein homolog B (S. cerevisiae)	
Immunogen Catalog Number: AG12981	Calculated MW: 25 kDa	

## Applications

Tested Applications:  
Cytometric bead array, Indirect ELISA

Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

# Selected Validation Data



Cytometric bead array standard curve of MP01774-1, SEC22B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85037-3-PBS. Detection antibody: 85037-4-PBS. Standard: Ag12981. Range: 0.313-40 ng/mL