

For Research Use Only

NME1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85038-3-PBS



Basic Information

Catalog Number: 85038-3-PBS	GenBank Accession Number: BC018994	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 4830	CloneNo.: 242670D1
Source: Rabbit	UNIPROT ID: P15531	
Isotype: IgG	Full Name: non-metastatic cells 1, protein (NM23A) expressed in	
Immunogen Catalog Number: AG1548	Calculated MW: 17 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01750-1, NME1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85038-3-PBS. Detection antibody: 85038-2-PBS. Standard: Ag1548. Range: 1.563-200 ng/mL.