

For Research Use Only

LAMB3 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85055-1-PBS



Basic Information

Catalog Number:

85055-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG25192

GenBank Accession Number:

BC075838

GeneID (NCBI):

3914

UNIPROT ID:

Q13751

Full Name:

laminin, beta 3

Calculated MW:

1172 aa, 130 kDa

Purification Method:

Protein A purification

CloneNo.:

242442F3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01812-2, LAMB3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85055-1-PBS. Detection antibody: 85055-3-PBS. Standard: Ag25192. Range: 0.781-100 ng/mL.