

For Research Use Only

# PTPN22 Recombinant antibody, PBS Only (Detector)

Catalog Number: 85062-1-PBS



## Basic Information

Catalog Number:

85062-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2397

GenBank Accession Number:

BC017785

GeneID (NCBI):

26191

UNIPROT ID:

Q9Y2R2

Full Name:

protein tyrosine phosphatase, non-receptor type 22 (lymphoid)

Calculated MW:

82 kDa, 92 kDa

Purification Method:

Protein A purification

CloneNo.:

242688A7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01818-2, PTPN22 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85062-2-PBS. Detection antibody: 85062-1-PBS. Standard: Ag2397. Range: 0.781-100 ng/mL.